QUIPSS II with Thin-slice TI₁ Periodic Saturation (Q2TIPS): A Method for Improving Accuracy of Quantitative Perfusion Imaging

W.-M. Luh, E. C. Wong*, P. A. Bandettini, and J. S. Hyde Biophysics Research Institute, Medical College of Wisconsin, Milwaukee, WI USA *Department of Radiology and Psychiatry, University of California, San Diego, CA USA

Introduction: QUIPSS II was introduced by Wong et al. (1) to improve quantification of perfusion using pulsed arterial spin labeling (ASL) techniques. A well-defined bolus of tagged blood is created and delivered to the imaging slices by applying a saturation pulse (90° pulse + crusher gradient) in the tagging region at TI₁ before the images are acquired at a later time TI₂ (Fig. 1a-b-c). Under these circumstances the perfusion signal is directly proportional to TI₁. In the previous implementation (2), residual errors occurred because of the imperfect match of the slice profile of the TI₁ saturation pulse to that of the inversion pulse at the distal edge (the edge further away from the heart), and the imperfect saturation of the TI₁ pulse. The present abstract describes and tests a method to minimize these errors. We refer to it as QUIPSS \underline{II} with $\underline{Thin-slice}$ $\underline{TI_1}$ Periodic \underline{S} aturation (Q2TIPS).

Methods: In Q2TIPS, the single saturation pulse at TI₁ of QUIPSS II is replaced by a series of smaller saturation pulses applied at the distal end of the tagging region. The distal edge of the saturation pulses match that of the inversion pulse. Fig. 1a-d-ef-c illustrates the timing and notation of this method where TI₁₅ is the stop time of the TI₁ saturation pulses. The TI₁ saturation slice is 2 cm and the TI₁ saturation pulse is a 3.2 ms sinc pulse applied every 50 ms and followed by a crusher gradient. All images were acquired using a 3T Bruker Biospec 30/60 scanner. A three-axis local head gradient coil and an endcapped quardrature birdcage coil were used. Single-shot blipped gradient echo EPI acquisition was used with 24 cm FOV and 8 mm slice thickness. The PICORE tagging scheme was used with a 10 cm tag created by a hyperbolic secant pulse with a 1 cm gap relative to the imaging slice. Tl₂ was 1.4 sec while TI₁ varied from 0.2-1.2 sec. TE was 27.2 ms and TR was 2.5 sec. Typically 100 repetitions were acquired and averaged after pairwise subtraction to generate a perfusion map. Perfusion quantification was performed using a kinetic model (3) with an assumed T₁ of blood of 1.5 sec.

Results: The thickness of the TI₁ saturation slice and the application rate of the TI₁ RF pulses determine the maximum flow velocity of the tagged blood to be saturated. The signal of the tagged blood with an average velocity less than 40 cm/sec traveling through the saturation slice will be destroyed. Since the tagged blood slows down as it approaches the capillary bed and the saturation slice is close to the imaging slice, current implementation should be sufficient to null the signal of remaining tagged blood which is destined for the imaging slices. When the velocity is higher than 40 cm/sec, it is more likely that blood will have passed the imaging slices. The minimum number of TI1 saturation pulses is determined by the time when the end of the tagged blood reaches the saturation slice. Although the TI, saturation pulses can be applied continuously until just before the imaging acquisition, excessive TI₁ saturation pulses will decrease the refreshing rate of the blood in the tagging region, which in turn increases the minimum TR. The anatomical image and the perfusion maps of QUIPSS II and Q2TIPS at TI₁ of 0.8 sec are shown in Fig. 2. The overall signal is lower in Q2TIPS and the bright signals in large veins also decrease. The averaged signals of QUIPSS II and Q2TIPS over gray matter ROIs are shown in Fig. 3a as a function of TI₁ and the calculated perfusion is shown in Fig. 3b with TI15 for Q2TIPS at 1.25 sec. The accuracy of perfusion measurement with Q2TIPS is dramatically improved as 15 evident from the linear regression curve with near-zero intercept in Fig. 3a. The calculated perfusion is overestimated in QUIPSS II mainly due to the imperfection of the TI₁ saturation pulse, especially at short TI₁. The averaged perfusion from the six data points obtained with Q2TIPS is 65±3 ml/100 ml/min. Fig. 3c shows the Q2TIPS signal as a function of TI_{1s} from 0.75 sec to 1.25 sec with TI₁ at 0.7 sec. For TI_{1s} too short, tagged blood at the tail of the tag reaches the imaging slice, and the inflection in this curve indicates the time at which the tail of the tag reaches the proximal end of the saturation band. Fig. 3d shows the Q2TIPS signal with TI₁/TI_{1s} at 0.7/1.1 sec as a function of TR. The signal decrease with TR shorter than 2.1 sec is due to the incomplete refreshing of blood in the tagging region. Since the TI₁ saturation slice is close to the imaging slice, the tagged blood destined for the imaging slice travels relatively slowly through the saturation slice and takes longer to replenish.

Conclusions: Q2TIPS improves the accuracy of quantitative perfusion measurement with pulsed ASL by providing a cleaner cut-off of the tag at TI₁, at the expense of slightly longer minimum TR.

References:

- 1. Wong, E. C., et al., ISMRM Proceedings, 85, 1997.
- 2. Wong, E. C., et al., ISMRM Proceedings, 1761, 1997.
- 3. Wong, E. C., et al., NMR Biomed. (in press).




Figure 1. Timing diagram of QUIPSS II (a-b-c) and Q2TIPS (a-d-e-f-c) as the tagged blood approaches the imaging slice. The diagram for the tagged blood represents only blood destined to go to the imaging slice.




Figure 2. (a) Anatomical and perfusion images of (b) QUIPSS II and (c) Q2TIPS with TI₁ at 0.8 sec.




Figure 3. (a) QUIPSS II and Q2TIPS signals and (b) calculated perfusion vz. TI₁. Q2TIPS signal vz. (c) TI₁₅ and (d) TR.

898 ADSTACL FORM TO SOSTING PRINCIPATIONS
INTERNATIONAL SOCIETY FOR
MAGNETIC RESONANCE IN MEDICINE
SIXTH SCIENTIFIC MEETING
SYDNEY AUSTRALIA
APRIL 18-24, 1998
ABSTRACT DEADLRE:
Received no liner than HOVEMBER 18, 1997

ABSTRACT DEADLINE:
Received no leter than NOYEMBER 18, 1997
All copyrights to accepted abstracts become the
property of the CURSUIT to property information
may be withheld by authors.

If the numbers have subminished than obstruct electromagnity in addition to their forms. Formula: MS WORD 5.00

Prefer Oral Presentation but willing to present as a power

Under required (as attable only for Oral Presentations).

Prefer Power but is illing to make Oral Presentation.

Protect outs

Protect outs

The authors certify that any work is thomas or animal subjects complete in the fire quality is due and of the relational or instrument and regulatory before Topose Categories to the one of the International or instrument and regulatory before Topose Categories to the one of Ill Into Doch it.

activate carrier resonant A (2) B (2) , C (1) , D (12) F (11)

Note the three and complete making ablesced the first asher Name Men-Ming Lub

Alberta Biophysics Res. Inst., Medical College of WT, B701 Watertoon Flank
Rd L Milwaukee, WI 53226

Commen U.S.A.

Telephone (414) 456-4592