



CEREBROVASCULAR CHANGES DURING THE VALSALVA MANEUVER MEASURED WITH VASO

4298




Daniel A Handwerker¹, Laurentius Huber¹, Puja Panwar¹, Benjamin Gutierrez¹, Javier Gonzalez-Castillo¹, Peter A Bandettini^{1,2}

1. Section on Functional Imaging Methods, LBC, National Institute of Mental Health, Bethesda, MD
2. Functional MRI Facility, National Institute of Mental Health, Bethesda, MD

INTRODUCTION

The Valsalva Maneuver involves an increase in chest pressure during a breath hold that alters blood pressure and heart rate (Dawson 1999). For a constant breath hold duration, the fMRI BOLD-weighted response scales with chest pressure magnitude. Both the undershoot during the breath hold and the signal peak after the hold ends increases with increasing pressure (Wu 2015). This parametric modulation of cerebrovascular reactivity by pressure has the potential to be used as a relatively simple cerebrovascular reactivity probe, but the mechanisms underlying the fMRI signal changes aren't fully understood. In this study, we use Vascular Space Occupancy (VASO) (Lu 2003) to measure BOLD-weighted and blood volume-weighted signal changes during Valsalva challenges with the goal of better understanding how the Valsalva Maneuver affects MRI signal changes.

VALSALVA TASK



VASO

SIEMENS Magnetom 7T, 32 channel Nova coil
Sequence: SS-SI VASO (Huber 2014): TR = 3 s, T11 = 1.2 s, T12 = 2.7 s
Inversion pulse: Tr-FOCI shape (Hurley et al., 2010), duration = 10 ms, bandwidth = 6.3 kHz
High-resolution protocol: 1.1x1.1x1.3 mm³, TE=24 ms, Partial Fourier=6/8
Low-resolution protocol: 2x2x2 mm³, TE=22ms, Partial Fourier=off
VASO might be contaminated by inflow of fresh (non-inverted) blood (Donahue et al., 2006). This is particularly challenging in hypercapnia related tasks because the arterial arrival time is reduced by up to 20% (Ho et al., 2011). These inflow effects are minimized here by the application of a custom designed inversion pulse with an adjustable inversion efficiency. (Here efficiency reduction = 7%.)
Readout: 3D EPI (Huber 2016b), GRAPPA = 2, Kernel size = 2 x 3, #reference lines = 24, ACS lines = FLASH (Talagala et al., 2016), phase correction = local,
Slab parameters: Slab excitation pulse duration 5 ms, variable flip angle 11°-90° to minimize across-segment smoothing, 10 slices/segments per slab.
To ensure proper spin inversion performance despite SAR and field inhomogeneity constraints at 7T, a MAFI B1+ map (Boulant, 2009) was acquired for every participant.

Data processing

Data were preprocessed using AFNI. Volumes were motion corrected within subject. The BOLD time series are the volumes without blood nulling. The VASO time series are the ratio of the nulled volume over the average of the flanking non-nulled volumes. Statistical significance maps (p<0.001 uncorrected) were separately calculated for BOLD and VASO data using a finite impulse response model. Since, SS-SI VASO data contain two different T1 weightings, a MP2RAGE-like Block model could be used to obtain quantitative T1-maps in functional EPI space (Huber, 2016a). Voxels with T1 values less than 4900ms that were significant in both BOLD and VASO were included in the clustering analyses. The average trial response for BOLD and VASO (18 volumes each) were separately scaled to mean 0 and standard deviation of 1. Ward's Hierarchical Clustering was run on these 36 volume vectors to identify clusters of voxels with similar BOLD and VASO response models.

Response Predictions

	Microvasculature		Large veins	
	BOLD	VASO	BOLD	VASO
[DeoxyHb]	↓	—	↓	—
CBF	↑	↓	↑	↓
CBV	↓	↓	↓	↓

An increase in blood pressure could potentially increase CBV and CSF volume

Since the relationship between deoxyhemoglobin concentration, blood flow (CBF), and blood volume (CBV) might vary, one cannot predict the VASO response from BOLD. For example, during a breath hold, there is an initial decrease in CBF and an increase in [DeoxyHb] which results in a BOLD decrease. In VASO, the larger vasoconstriction in microvasculature would cause a VASO increase while the increase in blood pressure might passively expand large veins and cause a VASO decrease


REFERENCES

Boulant et al., Counteracting radio frequency inhomogeneity in the human brain at 7 tesla using strongly modulating pulses, 2009, MRM, 61:1165-1172
Dawson et al., Critical closing pressure explains cerebral hemodynamics during the Valsalva maneuver, 1999, J Appl Physiol 86, 675
Donahue et al., Theoretical and experimental investigation of the VASO contrast mechanism, 2006, MRM, 56:1261-1273
Huber et al., Slab-selective, BOLD-corrected VASO at 7 tesla provides measures of cerebral blood volume reactivity with high signal-to-noise ratio, MRM, 2014, 72:137-148
Huber et al., Fast Dynamic Measurement of Functional T1 and Grey Matter Thickness Changes During Brain Activation at 7T, 2016a, ISMRM, p633.
Huber et al., Blood volume MRI with 3D-EPI-VASO: any benefits over SMS-VASO?, 2016b, p944
Hurley et al., Tailored RF pulse for magnetization inversion at ultrahigh field, 2010, MRM, 63:51-58
Lu et al., Functional magnetic resonance Imaging based on changes in vascular space occupancy, 2003, MRM
Talagala et al., Improvement of temporal signal-to-noise ratio of GRAPPA accelerated echo planar imaging using a FLASH based calibration scan, 2016, MRM
Wu et al., Effects of thoracic pressure changes on MRI signals in the brain, 2015, JCBFM


Supported by the NIMH Intramural Research Program

RESPONSE SHAPE PATTERNS


Sample Subject 2mm³ voxels (30 trials)




Sample Subject 1.15x1.15x1.7mm³ voxels (60 trials)



Breathing sensitive artifact in VASO



Motion artifact for SMS-EPI vs. 3D-EPI VASO



CONCLUSIONS

As seen in Wu et al 2015, the BOLD response usually decreases during the Valsalva breath hold and increases afterwards before returning to baseline
For voxels with similar BOLD responses, the VASO response varies more. For voxels in clusters with T1 values that are closer to gray matter, the VASO response also seems to decrease during the hold implying a gray-matter CBV increase.
With higher resolution voxels, we sometimes see a response spike in BOLD and VASO after the hold
Since this task has large, global blood flow changes and task-correlated chest movements, it has proven useful for highlighting and addressing artifacts that are hidden in more subtle studies.

ACKNOWLEDGEMENTS

We thank Benedikt Poser and Dimo Ivanov for sharing their nice 3D-EPI readout for the use in our VASO sequence. We thank Steve Cauley at MGH for sharing the interface of their image reconstruction for use with our SMS-VASO sequence as shared in the MGH CAIPI C2P
www.nmr.mgh.harvard.edu/software/c2p/sms
Paul Guillod wrote the stimulus presentation & pressure recording script in Presentation. Ronald Harper, George Dold and Thomas Talbot helped with the pressure recording system setup.

Portions of this study utilized the high-performance computational capabilities of the Biowulf Linux cluster at the NIH (http://biowulf.nih.gov)

Contact: dan.handwerker@nih.gov http://fim.nimh.nih.gov